

2020

TECHNOLOGY GUIDE



MY FACE } MY BODY
} 2019 {
AWARDS
DISTRIBUTOR
OF THE YEAR

MY FACE } MY BODY
} 2019 {
AWARDS
BEST CUSTOMER
SUPPORT

ABIA
Australian Beauty Industry Awards 2018
WINNER
2018 NATIONAL WHOLESALER
OF THE YEAR

ABIA
Australian Beauty Industry Awards 2019
WINNER
2019 NATIONAL WHOLESALER
OF THE YEAR



TheGlobalBeautyGroup
Technology for Beautiful Business

About Us

Welcome to The Global Beauty Group 2020 Technology Guide

This is a complete resource for Salon, Spa, Clinic and Practice owners looking to view our premier range of aesthetic technology. Here you will discover an exciting and forward-thinking selection of technology that greatly reflects our emphasis on high quality, innovation and excellent value.



OUR MISSION

what we do every day

We make people feel wonderful.

VISION OF SUCCESS

what we are striving to become

The Global Beauty Group is to be regarded by Salon, Spa, Clinic and Practice Owners as the No. 1 Aesthetic and Medical Grade Devices and Skin Care Company in the world.

OUR VALUES

what is non-negotiable

“Reliable, Accountable, Passionate, Innovative, Respectful, Supportive, Commercial and Fun”.



It's time to explore a new world of aesthetic technology

Simply select a topic from our menu below and flick through to its corresponding page. There you can view a great selection of devices designed to deliver performance-driven results for both your business and clients.

IPL Machines	4
Hair Removal Lasers	5
LED Devices	6
Microdermabrasion Machines	7
Oxygen Facial Technology	8-9
Microcurrent Machines	10
Radio Frequency Technology	11-12
Body Contouring Devices	13
Plasma Technology	14
Q-Switched Nd:YAG Lasers	15-16
Skin Analysis Systems	17
Intimate Rejuvenation Technology	18-19



IPL Machines

The Global Beauty Group's award winning range of multi-functional LUX Series IPL machines have been designed for results driven treatments, and are included on the Australian Register of Therapeutic Goods (ARTG No 295158).

SuperLUX and SuperLUX PRO

Award winning, next generation IPL technology

Winner of the 2017 Aestheticians Choice Award for the 'Multi-Function Equipment' category, SuperLUX SHR from LUX Series is a multifunctional, industry-leading technology.

- **Advanced IPL SHR Hair Removal** - The SuperLUX SHR handpiece offers rapid, pain free treatments for larger areas with a simple gliding motion that scans the treatment zone.
- **Superior and versatile results** - Integrated Dual-Mode technology combines IPL, radiofrequency (RF) and skin contact cooling to treat a wider range of skin types and hair colours, while a wide range of spot sizes enable rapid treatments on any body part.
- **Operator and client comfort** - Standard IPL treatments are fast and easy with both foot pedal and ergonomic handpiece firing options, plus -7° skin contact cooling keeps clients comfortable.
- **Powerful skin correction services** - Included filters allow clinicians to perform advanced Photo Rejuvenation, Acne Management, Pigmentation Correction and Vascular Therapy.
- **RF and Nd:YAG handpiece options** - Add laser tattoo removal and in demand skin tightening treatments with a LUX Series Radiofrequency (RF) skin tightening hand piece and Nd:YAG laser handpiece upgrade options.



Before



After 6 Treatments



Before



After 1 Treatment



UltraLUX and UltraLUX PRO

One device – 13 results driven treatments




Powerful and versatile, these multi-faceted devices allow beauty businesses to instantly add some of the most commonly requested treatments to their service menu:

- **Comprehensive range of IPL services** - With 3 IPL handpieces included, UltraLUX and UltraLUX PRO come complete with IPL SHR, IPL Dual-Mode for hair removal on lighter pigmented hair and your choice of spot size for standard IPL hair removal and IPL skin rejuvenation.
- **Ultrasound Body Contouring + RF Skin Tightening** - Slim and sculpt problem areas with Ultrasound Body Contouring, plus tighten lax or ageing skin on face + body with LUX Series RF.
- **Q-Switched Nd:YAG laser optional upgrade** - Add tattoo removal and carbon facials to your treatment menu to target new demographics with some of Australia's fastest growing treatments.

Before



After 5 Treatments



Hair Removal Lasers

Laser Hair Removal continues to provide a fast, easy and in demand solution for customers seeking effective permanent hair removal for virtually any part of their body.

HairLASE

The Smart Laser Hair Removal System

HairLASE is the latest innovation to join the exclusive LUX Series range of treatment technologies:

- **Rapid, hassle free treatments** - HairLASE offers default STAMP and MOTION preloaded settings, plus easily adjustable pulse width options for fast and effective treatments on Fitzpatrick Skin Types I-VI.
- **Superior patient comfort** - Advanced skin contact cooling of -4°C sets a benchmark for exceptional patient comfort.
- **Streamlined design** - A smart compact system that integrates an ergonomic handpiece and the best of diode laser technology, HairLASE fits seamlessly into all treatment rooms.
- **Pre-programmed treatments** - With a modern, user-friendly LCD touch screen interface, HairLASE recommends best settings for each client skin type.



LED Devices

With its restorative properties, hands-free application and visible results, LED has become a powerful clinical tool for promoting overall skin and cellular wellbeing with little time, effort and consumables required.

MediLUX

Ultimate versatility for unprecedented LED results

Combining 5 LED wavelengths, wide adjustable head for maximum coverage and one of the highest power outputs available, MediLUX allows clinicians to push their LED service potential further than ever before:

- **Pre-programmed LED treatments** - Easy to use touch-and-go facial services on a 10.4 inch colour touchscreen interface.
- **Manual mode treatments** - Ultimate versatility with customisable time, sequence and power options on blue, green, yellow, red and near-infrared light settings.
- **Adjustable LED panel** - Panel adjusts to the curves of the face and neck or to completely straight for unique, non-facial services such as the back and stomach.
- **Near-Infrared services** - Boost adenosine triphosphate production (ATP) to promote wound healing, release of toxins, DNA repair and cellular regeneration.
- **Ultimate power output** - Higher wavelengths for faster treatments in Acne Management, Pigmentation Correction, Redness Reduction, Age Management and Dermal Repair.



Bio-Synthesis™

Industry-leading multi wavelength light therapy

Bio-Synthesis™ adds the professional edge to any skin corrective service:

- **4 wavelengths of light therapies** - Features blue light for skin clearing, red light for age management, yellow light for calming inflammation and green light for skin tone correction.
- **Pre-programmed and manual modes** - Tailor treatments to your clients' needs with adjustable parameters.

Before



After 2 Treatments



bt-Accent LED™ and Bio-Penta™

The sleek, compact and affordable LED solution



This multiple award winning device equips both beginner and advanced aestheticians with straightforward, results-driven light activated facials:

- **3 core LED therapies** - Features blue light for skin clearing, red light for age management, and purple light for dual-action skin clearing, age management and skin balancing services.
- **Easy pre-set treatment times** - Select from 10, 20 or 30 minute treatment times to provide comprehensive LED facials.
- **Lightweight, compact design** - Streamlined construction allows easy mobility between treatment rooms, while still providing excellent treatment coverage to the face, neck and décolleté.
- **Hydro/microdermabrasion, ultrasonic exfoliation and product infusion options** - Combining the advanced hydrodermabrasion and microdermabrasion of Bio-Hydroderm Trinity™ with the blue, red and combination light therapy of bt-Accent™ LED, Bio-Penta™ also utilises microsonic facial cleansing and microcurrent + ultrasound product infusion for the ultimate 5-in-1 skin care system.

Micro/Hydrodermabrasion Machines

The ability to provide clients with smooth, thoroughly exfoliated, hydrated and youthful looking skin in as little as one treatment is a key factor in the ongoing popularity of microdermabrasion and hydrodermabrasion. Using clinical strength techniques to exfoliate and renew the skin, treatments can also be combined with skin corrective infusion products for results that are immediately psychologically satisfying and promote long-term skin health.

HydroLUX

Advanced Skin Resurfacing + Product Infusion



Combining the best of both microdermabrasion and hydrodermabrasion, HydroLUX delivers unprecedented service versatility and treatment scope.

- **Customised treatment options** - Select from 10 levels of vacuum and flow across 3 easily interchangeable 300ml solution containers to tailor treatments to your client's unique skin concerns.
- **Interchangeable treatment tips** - Choose from an 8mm, 15mm and 25mm tip size in an 80, 100, 150, 200 and 250 grit intensity for customised treatments.
- **One touch cleaning** - Cleaning is easy, fast and effortless with an inbuilt one press 2-minute cleaning cycle.
- **Wet to dry microdermabrasion** - Switch from dry microdermabrasion to liquid infusion hydrodermabrasion with a press of the footswitch, making dry passes on problem areas easier than ever before.

Oxygen Facial Technology


The next generation of oxygen facial technology combines oxygenating for a radiant glow with the infusion of innovative ingredients for enhanced skin corrective results. This also prepares the skin for further technology enhanced treatments that work in synergy with oxygenation to create instantly refreshed, rejuvenated and younger looking skin.

Geneo™

The ultimate 3-in-1 super facial

Oxygenate and invigorate tired, stressed and congested skin with the industry's most innovative oxygen therapy technology:

- **Patented Geneo™ technology** - Geneo™ application generates a chemical reaction producing CO2 bubbles which gently burst on the skin's surface to revitalise skin.
- **Simultaneous exfoliation and ultrasound infusion** - Anti-ageing formulas are infused into the skin via oxygen and ultrasound while the unique Oxypod tip exfoliates for a smoother, fresher and plumper complexion.
- **Radiofrequency (RF) skin tightening and product infusion available** - TriPollar® massager handpiece delivers additional skin tightening services to the face, neck and hands to further enhance the anti-ageing effects.



Geneo™ Oxypods

Innovative and Unique Ingredients

Specially developed Oxypods for Geneo combine advanced research backed ingredients like tree bark extract, amino acid derivatives, salicylic acid, jojoba, argan oils, edelweiss, vanilla and horsetail extracts. When combined with Geneo Primer Gel during treatment, the uniquely formulated products provide nutrition to the skin via a synergistic combination of active ingredients for 3 main benefits: exfoliation, oxygenation and nourishment.

Geneo™ Illuminate

Vitamin C Facial

- Geneo Illuminate unifies skin tone, improves pigmentation and rejuvenates the skin.
- The Geneo Illuminate OxyPod includes Kojic Acid and Vitamin C for an even skin tone effect, as well as Opuntia fruit extract and Carrot Extract for moisturizing and smoothening the skin.
- The Geneo Illuminating Serum completes the treatment effect with powerful brightening ingredients such as Tree Bark Extract and Amino acid derivatives, and moisturizing ingredients and nutrients like Triglycerides, Jojoba and Argan oil complex, and Vanilla and Horsetail extract complex.

Geneo™ Revive

Red Algae Facial

- Geneo Revive refreshes dull skin, reduces fine lines and improves overall skin texture and appearance.
- The Geneo Revive OxyPod includes Red Algae a powerful antioxidant with an anti-aging effect, and also Caffeic Acid and Retinol, known to treat signs of aging and assist in skin rejuvenation.
- The Geneo Reviving Serum completes the treatment effect with additional ingredients that encourage collagen production and reduce signs of aging, like Peptides complex and Edelweiss, as well as nutrients and vitamins such as Jojoba and Argan oil complex, Vanilla and Horsetail extract complex.

Geneo™ Balance

Charcoal Facial

- Geneo Balance treatment purifies and soothes oily skin and improves its texture and appearance.
- The Geneo Balance OxyPod includes Bamboo Charcoal for a detoxifying effect, Camu-Camu extract for vitamins and antioxidants, and Mandelic Acid for treating acne-prone skin.
- The Geneo Balancing Serum completes the treatment effect with powerful antioxidants, vitamins, nutrients and purifying ingredients such as Salicylic acid, Peptides, Jojoba and Argan oils, Vanilla and Horsetail extracts and Potassium Azeloyl Diglycinate.

Microcurrent Machines

Lifted, firmed and contoured facial features can now be achieved without invasive procedures thanks to microcurrent technology. Also known as “non-surgical facelifts” clients of all skin types and ages can enjoy a younger, fresher appearance without pain, injectables or downtime.

bt-Nano™

Electrify your anti-ageing treatments



Compact size, maximum efficacy — bt-Nano™ is the non-surgical facelift you can take anywhere:

- **5 pre-set facial services** - Select from 5 facial service modes including the bt-Cocktail™ Lift, Quickly Clear Facial™, Platinum Touch Hand Renewal™ and more.
- **Lightweight, portable design** - Battery operated and weighing less than 500 grams, this compact device is perfect for spas and clinics with multiple treatment rooms, and is also ideal for aestheticians on the go.

Bio-Ultimate Platinum™

The industry's leading non-surgical facelift



Lift, tone, firm, contour and volumise with this high-powered, sophisticated and feature-packed microcurrent system:

- **15+ pre-programmed facial services** - Deliver ultra-professional facial contouring services tailored to your client's specific skin concerns and treatment objectives, from age defying treatments to skin clearing, turbo lifting facials, and more.
- **Patented Suzuki Sequencing™** - Premium engineered system enhanced by patented Suzuki Sequencing™, and complete with full facial and hand accessories for comprehensive anti-ageing treatments.
- **Stand-alone or layered services** - Bio-Ultimate Platinum™ delivers powerful stand-alone facial lifting and contouring services, but can also be used as part of a layered technology treatment plan, such as post-hydrodermabrasion/microdermabrasion.

Before



After 6 Treatments



Radiofrequency Technology

Delivering a wide range of treatment applications, radiofrequency technology (also known as RF) generates collisions between charged molecules and ions to create heat energy. RF energy is commonly used in aesthetic treatments to produce a microinflammatory stimulation of fibroblasts which helps to generate collagen and elastin, as well as to effectively stimulate fat breakdown, induce cellulite reduction, and tighten lax skin.

Pollogen MAXIMUS™

Clinically proven body contouring

Sculpt, contour and rejuvenate with the complete anti-ageing and body contouring platform, Pollogen MAXIMUS™:

- **TriLipo® with DMA™ (Dynamic Muscle Activation) technology** - Reduce circumferential fat and cellulite by heating and pressurising the targeted subcutaneous fat to the point of release. Included DMA™ technology provides simultaneous muscle lifting and toning for more dramatic body contouring and figure sculpting results.
- **Multiple treatment benefits** - TriPollar® clients can also enjoy dermal volumising, skin tightening, stretch mark reduction, wrinkle reduction, cellulite reduction and body detoxification.
- **3 treatment applicators** - A variety of applicator sizes are provided to deliver effective accessibility and treatment coverage to delicate, small-medium, and large body areas.
- **No discomfort or downtime** - Patients experience comfortable, convenient in-and-out treatments without disruption to work or daily activities.



Before

After 10 Treatments



Before

After 5 Treatments



Divine Pro™

Multi-dimensional facial rejuvenation



The Divine Pro™ is a multi-solution platform that works to non-surgically resurface and rejuvenate both the epidermal and dermal layers of skin:

- **VoluDerm™ RF enhanced microneedling** - Using minimally invasive micro needles (0.5mm H x 150 micron W), the VoluDerm™ handpiece engages Radiofrequency (RF) technology to renew the epidermal layer and volumise the dermal layer.
- **TriFractional™ for epidermal rejuvenation** - Using minimally invasive pins (0.2mm H x 300 micron W), TriFractional™ treats the outer epidermis, functioning as a mild peeling procedure to improve skin tone, texture and clarity.
- **TriPollar® with DMA™ technology** - For high-tech anti-ageing facials, TriPollar® RF technology enhanced with DMA™ techniques tighten, lift, contour, volumise and induce collagen remodelling to the face, neck and décolleté.
- **Cutting-edge innovation** - Divine Pro™ delivers a holistic approach to anti-ageing by repairing the skin from the inside out.
- **Minimal discomfort and downtime** - Treatments range from non-invasive to minimally invasive, meaning clients experience superior comfort and decreased downtime when compared to other skin resurfacing technologies.
- **Quick 30 minute sessions** - Treatment takes on average 30 minutes per session, offering absolute convenience for clients.

Before



After 5 Treatments



Before



After 4 Treatments



Body Contouring Devices

With 67% of Australians now overweight or obese (Australian Institute of Health and Welfare, 2019) body sculpting treatments are now one of the fastest growing aesthetic treatments amongst both men and women.

Criocuum

Professional Cryo Body Contouring for results-driven therapists

Experience the latest innovation in Cryo Body Contouring technology:

- **Advanced Cryo Body Contouring technology** - Adjustable temperatures as low as -10° are applied to the treatment area to slim and sculpt, with visible results as early as 15 days following treatment.
- **Versatile treatment areas** - With large fold, small fold, mini fold and curved contour handpieces clinicians can treat a wide range of areas including the chin, upper arms, abdomen, thighs, buttocks, back and hips using their choice of 1 or 2 handpieces at a time.

Before



After 3 Treatments



LUX Series

The pain-free approach to body enhancement

Slim, tone, tighten and contour with SlimLUX, UltraLUX and UltraLUX PRO devices:

- **Results driven services** - A non-surgical and non-invasive alternative to liposuction and face/skin lifting procedures with visible results after as little as 1 treatment.
- **No downtime** - Lifestyle-driven benefits include quick, comfortable treatments, and no post-treatment downtime, so clients are free to resume daily activities as normal.

Before



After 2 Treatments



Plasma Technology

An increasingly popular non-surgical cosmetic dermatology technique, plasma treatments utilise a plasma flash (small electrical charge) to create a radius of microplasma that delivers energy to the superficial skin cells. This transitions the targeted treatment area from a solid to a gas state without damaging surrounding tissue (a process known as sublimation to initiate a wound healing response in the skin, stimulating fibroblast cells to produce more collagen and repair the skin.

Plexr Plus

The next generation of medical-grade plasma technology



Non-surgical procedures to smooth, tighten, lift, repair and rejuvenate the face and body:

- **Advanced skin rejuvenation results** - Treatments deliver results on frequent client face and body concerns; including fine lines and wrinkles, loose and lax skin, drooping eyelids and epidermal scarring.
- **Innovative treatment technology** - A 500 micron treatment tip delivers energy only to superficial skin cells, without involving the basal layer. This reduces risk and side effects associated with more invasive skin rejuvenation treatments.
- **Non-surgical blepharoplasty** - Plexr Plus is a clinically proven treatment for lifting drooping skin on the upper and lower eyelids, resulting in a reduction of excess skin in the external region around the eye and a tightening of lax skin and wrinkles.
- **Skin lesion removal** - Clinicians can provide effective removal of moles, skin tags, senile warts and cherry angiomas on almost all areas of the face and body.



Q-Switched Nd:YAG Lasers

By producing short, nanosecond long pulses of powerful energy, Q-Switched Nd:YAG lasers are able to safely treat the particles within skin pigmentation and tattoo inks. With 34% of tattooed Australians either regretting a tattoo or considering having it removed (Choice, 2017), tattoo removal services are a popular treatment for salons, spas, clinics or practices wanting to expand their client base.

UltraLASE

Laser technology for unwanted tattoos and skin rejuvenation

This dedicated Q-Switched Nd:YAG laser system from LUX Series allows clinicians to remove unwanted tattoos and perform the highly sought after Carbon Facial, also known as the “China Doll Facial” or “Hollywood Facial”:

- **Fast and easy treatments** - Capable of removing tattoos containing black, dark blue, green and red pigments within 6–12 sessions depending on the age and quality of the tattoo.
- **On trend Carbon Facials** - Carbon creme is applied to the entire facial area which is absorbed into the pores. This is followed by application of laser light wavelengths to activate the carbon, shattering it into particles to clean and rejuvenate the epidermal layer of the skin.
- **Excellent clarity and ease of use** - Ultra-bright aiming beam offers maximum clarity and ease of application while the easy to use interface enables seamless switching between 532nm and 1064nm wavelengths.
- **Multifunctional options available** - The UltraLASE Q-Switched Nd:YAG laser handpiece is also an optional upgrade on the UltraLUX PRO and SuperLUX PRO multifunctional systems.



Before



After 4 Treatments



After 6 Treatments



Before



After 4 Treatments



After 7 Treatments



TRI-BEAM™

Industry-leading laser technology to erase, rejuvenate, renew and regenerate



Included on the ARTG (No. 223320), TRI-BEAM™ is the world's fastest tattoo removal laser. It's advanced Q-Switched Nd:YAG laser boasts outstanding skin resurfacing and tattoo removal abilities:

- **Full colour spectrum** - 532nm, 585nm, 650nm and 1064nm wavelengths allow for comprehensive treatment of inks across the colour spectrum at up to 20hz pulse capacity.
- **Flexible spot sizes and rapid delivery** - Adjustable spot sizes (2mm to 10mm) facilitate maximum coverage.
- **Broad range of applications** - A 532nm wavelength treats freckles, solar lentigines and superficial pigmentation.
- **Deep laser penetration** - A 1064nm wavelength targets post inflammatory hyper-pigmentation, melasma, nevus of ota and ABNOM by shattering the pigment, without thermal trauma.
- **Exceptional Carbon Facial results** - A 1064nm Gen Technique™ laser promotes collagen regeneration to thermally lift, tighten and tone the skin, plus reduce the appearance of fine lines.
- **Reduce scarring and regenerate skin** - Promote a healthier, more youthful complexion by stimulating the production of collagen and elastin to improve overall skin tone, reduce scarring and decrease pore size.
- **Thorough energy distribution** - The advanced beam profile ensures an even delivery of energy to tissue throughout the pulse, making for unrivalled light tissue interaction.

Before



After 1 Treatments



Before



After 6 Treatments



Skin Analysis Systems


Skin analysis forms an integral part of any client consult. Advanced analysis technology can assist with accurately diagnosing surface and subsurface skin concerns, enabling aestheticians to educate clients about their skin and confidently proceed with the right skincare or treatment for their needs.

OBSERV 520

The first step to optimum skin health

Take your clients on a revolutionary journey beneath the skin's surface:

- **State of the art skin observation** - OBSERV 520 assists skin professionals with accurate skin condition diagnosis at both the epidermal and dermal layer using patent-pending skin fluorescence technology. This scientific principle exposes the skin to low, safe levels of UV light, highlighting skin characteristics and patterns with unparalleled clarity when compared to similar diagnostic tools.
- **6 observation modes** - Clinically examine the entire facial region with Daylight mode, True UV mode, Simulated Wood's lamp mode, Parallel-polarisation mode, Cross-polarisation mode and Complexion analysis.
- **Impressive side by side comparisons** - Keep clients engaged on their treatment journey with cutting-edge side by side comparisons that track improvements and progress.
- **Modern, user-friendly system** - Control the OBSERV 520 system conveniently and wirelessly using an iPad and the downloadable OBSERV 520 app to track clients skin journeys.
- **Advanced Imaging System (AIS)** - A rotating chinrest allows clinicians to take full-face digital images front on, and at angles of 40 and 90 degrees. FaceNotes feature allows clinicians to write and draw directly on examination images to highlight areas of concern, then print or email them directly to the client with their skin analysis report.
- **Customised to your business** - Easily access your treatment and skin care options within the consultation, for an organic transition from skin examination to treatment recommendations and solid bookings.



Intimate Rejuvenation Technology

In a study of 8,081 postmenopausal women, 38% experienced symptoms of vaginal atrophy including urinary incontinence, vaginal inflammation and discomfort during intercourse (REVIVE Survey, 2012). With advances in treatment technology, intimate rejuvenation gynecological treatments are becoming increasingly popular in clinics that support their clients achieve a happier and healthier quality of life.

Thuzzle

Combined thermal energies that solve the “puzzle of beauty”

Complete rejuvenation of the face and body with one innovative, all-in-one device:


- **In demand therapeutic treatments** - Comfortable intimate rejuvenation treatments utilise monopolar and bipolar modes to regenerate collagen. This tightens the vaginal area elasticity to address commonly experienced concerns such as incontinence, sexual discomfort and weakening of the pelvic floor.
- **Aesthetic applications** - Radio Frequency (RF) can be utilised for Vulva Tightening, to minimise the labia minora and labia majora to increase client confidence in the appearance of intimate areas.
- **A simple, hygienic solution** - Disposable vaginal handpiece is designed for 1-person use and is suitable for an entire patient treatment series (6-8 treatments).
- **Safe and comfortable** - From industry leader GMV, Thuzzle includes inbuilt handpiece temperature control, motion control and skin contact recognition for treatments that are easy, comfortable and safe.
- **Client results** - Studies show a 30% performance improvement on patient outcomes compared to treatments with other thermal energy devices available.




Before



After 5 Treatments



Before



After 4 Treatments



Before



After 2 Treatments



Before



After 1 Treatment



- **Exclusive technology** - Thuzze is the only fractional technology capable of delivering a radial type signal to treat tissue uniformly, and is uniquely capable of delivering up to four electromagnetic waves controlled at four different depths for enhanced treatment outcomes.
- **Complete clinic solution** - From fractional facial rejuvenation technology to control crosslinked double, triple and tetra wave handpieces, Thuzze is the ideal solution for full face and body rejuvenation. This all-in-one device provides some of the industry's fastest growing treatments including advanced body contouring, scar reduction, skin tightening, skin rejuvenation, intimate rejuvenation and lymphatic drainage.



50%

Increase in
non-surgical
fat reduction



29%

Increase in
micro-ablative
fractional resurfacing



53%

Increase in
vaginal
procedures



12%

Increase in
non-surgical facial
rejuvenation services

** 2018 report by The American Society for Aesthetic Plastic Surgery (ASAPS) into aesthetic treatments identified demand increases between 2014-2018*

Why do 2,200+ salons, spas, clinics and practices across Australia and New Zealand choose The Global Beauty Group?

"The support we have received from The Global Beauty Group has been fabulous. We have purchased multiple technologies through your company now and have been really happy with all of them. The initial training and ongoing support has also been great!"

– Rikki Ballantyne, Skin Forum

"Time and time again, I keep coming back to The Global Beauty Group because they always have the best technology to fit the purpose and our needs"

– Kerree Kerr, Self-Centred Medi Spa

Our IPL machine is a real workhorse - we get fantastic results. Training is awesome! Follow up is always punctual. When it's time to upgrade or upscale I will give the team at GBG a call again"

– Gry Tomte, HÜD Skin + Body

"Partnering with The Global Beauty Group has been really easy. It was seamless. The one on one training & marketing help has been incredible!"

– Olivia Doncaster, Olivia Catherine Skin Concept

"Reliable machines and excellent communication. Would highly recommend"

– Anita Turner, Gippsland Cosmetic Laser Clinic

"Our technology allows superior long-term treatment and retention of patients. It allows better and long term engagement with patients. For me personally, it has allowed us to slowly bring patient's conditions under control without the downtime"

– Associate Professor Greg Goodman, Dermatology Institute of Victoria



PROUD SUPPLIERS TO

RESULTS LASER CLINIC

BRAZILIAN BEAUTY

dermalogica

Brazilian Butterfly

NATIONAL ACADEMY OF BEAUTY

DERMATOLOGY INSTITUTE OF VICTORIA

Clearskincare Clinics

endota spa

cosmetic

WEST COAST INSTITUTE

Ella Baché

SILK LASER CLINICS

DANIEL LANZER

Australian Skin Clinics

NSW TAFE